



## **PURELAB® Classic UF\***





## DI Pack

ELC186

Lifetime expectancy\*\*: 12 months What is purified from the water? lones Impact of not changing: Weakly-bound ions released back into the water



Ultra Micron Filter ELC151 Lifetime expectancy\*\*: 24 months What is purified from the water? Bacteria and microbial by-products such as endotoxins Impact of not changing: Bacterial and endotoxin specification will not be maintained reduced flow rate



Point Of Use Filter ELC134 Lifetime expectancy\*\*: 6 months What is purified from the water? Particulates and bacteria Impact of not changing? Potential recontamination of purified water and reduced flow rates

\* A combination of treatments and technologies produce the final water specification

\*\* Service life will depend on the application and feed water quality

