



PURELAB® Classic UV*





DI Pack

ELC186

Lifetime expectancy**: 12 months What is purified from the water? lones

Impact of not changing:

Weakly-bound ions released back into the water



UV Lamp (185 nm)

ELC118

Lifetime expectancy**: 12 months What is purified from the water? Purifies TOC, inactivates DNA & RNA Impact of not changing:

Uncontrolled bacterial contamination TOC levels not reduced



Point Of Use Filter***

ELC134

Lifetime expectancy**: 6 months What is purified from the water? Particulates and bacteria Impact of not changing?

Potential recontamination of purified water and reduced flow rates

- * A combination of treatments and technologies produce the final water specification
- ** Service life will depend on the application and feed water quality
- *** For DNA, RNA & endotoxin specification at point of use, please use a ELC197 Biofilter

