



PURELAB® flex 1*





DI Pack

ELC208

Lifetime expectancy**: 12 months What is purified from the water? lons Impact of not changing: Weakly bound ions are released back into the water reducing resistivity



Point Of Use Filter*** ELC145

Lifetime expectancy**: 3 months What is purified from the water? Particulates and bacteria Impact of not changing: Potential recontamination of purified water and reduced flow rates

* A combination of treatments and technologies produce the final water specification

** Service life will depend on the application and feed water quality

*** For DNA, RNA & endotoxin specification at point of use, please use a ELC197 $\,$

