



# PURELAB® flex 2\*



## DI Pack

ELC208

Lifetime expectancy\*\*: 12 months

What is purified from the water?

Ions

Impact of not changing:

Weakly bound ions are released back into the water reducing resistivity



## UV Lamp 185/254 nm

ELC210

Lifetime expectancy\*\*: 12-18 months

What is purified from the water?

Purifies TOC, inactivates DNA and RNA

Impact of not changing?

Uncontrolled bacterial contamination and TOC levels are not reduced



## Point Of Use Filter\*\*\*

ELC145

Lifetime expectancy\*\*: 3 months

What is purified from the water?

Particulates and bacteria

Impact of not changing:

Potential recontamination of purified water and reduced flow rates

\* A combination of treatments and technologies produce the final water specification

\*\* Service life will depend on the application and feed water quality

\*\*\* For DNA, RNA & endotoxin specification at point of use, please use a ELC197

