



Think inside
the box

PURELAB® Quest



DI Pack ELC292

Lifetime expectancy*: 6 months
What impurities are removed from the water?
 Ions
Impact of not changing:
 Weakly bound ions are released back into the water reducing resistivity.



Reverse Osmosis ELC217

Lifetime expectancy*: 24 months
What impurities are removed from the water?
 Majority of organics and inorganics
Impact of not changing: High levels of organics and inorganics remain in the water reducing the functional life of the DI Pack.
 Note: Replacement need technician support.



Vent Filter ELC291

Lifetime expectancy*: 12 months
What impurities are removed from the water?
 Protects the sealed system against bacteria/particulates in the air
Impact of not changing:
 The filter can potentially become blocked reducing its function to adequately allow air through.



Point Of Use Filter** ELC134

Lifetime expectancy*: 3 months
What impurities are removed from the water?
 Particulates and bacteria
Impact of not changing:
 Potential recontamination of purified water and reduced flow rates.

* Service life will depend on the application and feed water quality

** For DNA, RNA & endotoxin specification at point of use, please use a ELC197 Biofilter

