



MEDICA PRO LPS*





DI Pack Lifetime expectancy**: 3-12 months

ELC174

What is purified from the water? lons Impact of not changing: Weakly bound ions are released back into the water reducing resistivity



Composite Vent Filter ELC136

Lifetime expectancy**: 6 months What is purified from the water? Reduces CO₂ and airborne organics from entering the sealed system Impact of not changing: Increase in airborne contaminants and potential recontamination of purified water



Pre-Treatment L3

Lifetime expectancy**: 6 months What is purified from the water? lons

Impact of not changing:

Chlorine, chloramine and organics remain increasing the burden on downstream consumables which will need more frequent and expensive changes



UVC Lamp 254 nm

ELC105 Lifetime expectancy**: 12 months What is purified from the water? Inactivates DNA and RNA

Impact of not changing: Uncontrolled bacterial contamination

ELC109

ELC192



Ultra-Micron Filter Lifetime expectancy**: 12 months

What is purified from the water? Particles & microorganisms Impact of not changing: Reduced flow-rates and bacterial specification will not be maintained

* A combination of treatments and technologies produce the final water specification

** Service life will depend on the application and feed water quality

