



Think inside
the box

PURELAB® Chorus 3*



Desinfection Tablets ECT1

Lifetime expectancy:** 6 months
What is purified from the water?
Sanitizes tubing or pipe work to remove bacteria, reduce bio film and protect the system against bacterial growth.



Pre-Treatment ELC241

Lifetime expectancy:** 6 months
What is purified from the water?
Chlorine and particulates
Impact of not changing? Chlorine, chloramine and organics remain increasing the burden on downstream consumables which will need more frequent and expensive changes



Reverse Osmosis ELC240

Lifetime expectancy:** 24 months
What is purified from the water?
Majority of organics and inorganics
Impact of not changing?
High levels of organics and inorganics remain in the water reducing the functional life of the DI Pack (Pre-Treatment)
Note: Replacement need technician support



Composite Vent Filter ELC216

Lifetime expectancy:** 6 months
What protection does it provide?
Reduces CO₂ and airborne organics from entering the sealed system
Impact of not changing?
Increase in airborne contaminants and potential recontamination of purified water

* A combination of treatments and technologies produce the final water specification

** Service life will depend on the application and feed water quality

